



## Pulse oximeter PM 100 connect

Art.-No. 79456 | EAN 4015588 79456 8

- Bluetooth® data transfer to the VitaDock+ app for iOS and Android and to VitaDock Online\*
- Pulse oximeter to measure the saturation of oxygen in blood (SpO2 %) and the heart rate (pulse)
- View measurements on a OLED display
- Simple „one touch“ operation – ideal for on the go



Smartphone not included



### Application | Features

- Bluetooth® data transfer to the VitaDock+ app for iOS and Android and to VitaDock Online\*
- Pulse oximeter to measure the saturation of oxygen in blood (SpO2%) and the heart rate (pulse)
- View measurements on OLED display
- Simple „one touch“ operation – ideal for on the go
- Oxygen saturation (SpO2) and pulse rate can be displayed as bars or waves
- Automatic switch-off after 8 seconds
- Low battery display
- Incl. carrying strap
- Incl. batteries (2 × 1.5V AAA)
- Certified medical device\*\*

### Technical Data

- Size: approx. 5.8 × 3.4 × 3.5 cm
- Weight: approx. 0.055 kg
- Measuring range: SpO2: 70% ~ 100%
- Pulse: 30-235 beats/min.
- Batteries: 2 × 1.5V (AAA)

### Scope of Delivery

- Pulse oximeter
- Batteries
- Instruction manual

### Logistical Data

- Sales packaging**  
1 pieces | 12.6 × 18 × 6.6 cm | 0.3 kg
- Packaging unit**  
4 pieces | 37 × 14.5 × 13 cm | 1.2 kg

\*Bluetooth® 4.0 compatible with VitaDock+ App as of Apple (iOS 7.0), iPhone® 4s, iPod Touch® (as of 4 Gen.), iPad 3, as of Android 4.3  
\*\* Certification according to MDD (Medical Device Directive)



06\_23SEP2021