



- ✓ Used to help relieve bruises, swelling, muscle spasms, aches and pains, headaches and first-degree burns
- ✓ Ultra-soft, flexible fabric for optimal fit and comfort
- ✓ Durable and reusable



Hot & Cold Pack GP 101

Reusable gel compress

Art. No. 88101 | EAN 40 15588 88101 5

Filled with soft,
form-fitting hydrogel

APPLICATION | FEATURES

- Treatment area: neck
- For hot and cold usage
- Cold therapy: place the product in freezer or refrigerator for at least 2 hours
- Hot therapy: heat the product in microwave for 35 seconds (800 watts)
- Apply hot/cold product for no more than 15 minutes at a time
- Remains flexible even when frozen for comfortable use (-18 degrees)
- Easy to clean fabric
- Non-toxic
- Expiration date: 5 years
- Certified medical device (MDR Class 1)

TECHNICAL DATA

- Size: approx. 52 x 18 cm
- Weight: approx. 494 g
- For a neck circumference from 34 to 52 cm

SCOPE OF DELIVERY

- Hot & Cold Pack
- Instruction manual

LOGISTICAL DATA

Sales packaging
1 pc | 9.2 x 9.7 x 20.3 cm | 0.7 kg

Packaging unit
4 pcs | 37 x 19 x 11 cm | 2.8 kg

03AUG2026